


MASSTABELLE

mediven elegance[®], mediven for men[®], mediven comfort[®],
 mediven plus[®], duomed basic[®]

Standard ohne Naht							
Umfänge	I	II	III	IV	V	VI	VII
cG Oberschenkelstrumpf mit Haftband	43 – 48	45 – 52	49 – 56	53 – 60	56 – 64	60 – 68	64 – 72
cG Oberschenkelstrumpf mit Haftband extra weit	49 – 57	53 – 62	57 – 67	61 – 72	65 – 77	69 – 81	73 – 85

	cG	43 – 57	45 – 62	49 – 67	53 – 72	56 – 77	60 – 81	64 – 85
	cF	39 – 52	41 – 56	44 – 60	47 – 65	50 – 69	53 – 73	57 – 78
	cE	30 – 37	33 – 40	35 – 43	37 – 45	39 – 48	41 – 51	42 – 53
	cD	27 – 33	29 – 36	32 – 39	34 – 42	36 – 45	38 – 48	40 – 50
	cC	28 – 34	30 – 37	33 – 40	35 – 43	37 – 46	39 – 49	41 – 51
	cB ₁	23 – 27	24 – 29	26 – 32	29 – 35	31 – 37	33 – 39	35 – 41
	cB	18 – 20	20 – 22	22 – 24	24 – 26	26 – 28	28 – 30	30 – 32
	cY	26 – 31	28 – 33	29 – 35	31 – 37	32 – 38	33 – 40	34 – 42
	cA	17 – 22	19 – 24	21 – 26	23 – 29	25 – 32	27 – 34	28 – 35

Länge	
/G normal	72 – 83
/G kurz	62 – 71
/D normal	39 – 44
/D kurz	34 – 38

Maxi-Leibteil AT mediven plus		
Grösse I – III	T = ab 95 cm	H = ab 110 cm
Grösse IV – V	T = ab 120 cm	H = ab 135 cm
Grösse VI – VII	T = ab 140 cm	H = ab 145 cm


Cosanum AG
 Brandstrasse 28
 8952 Schlieren
 Kundenservice Fachhandel/ Arzt:
 Tel. 043 433 66 77

Gratis Bestellfax:
0800 730 900

MASSTABELLE

mediven forte®

Standard ohne Naht							
Umfänge	I	II	III	IV	V	VI	VII
cG Oberschenkelstrumpf mit Haftband	44 – 47	46 – 50	50 – 54	53 – 58	56 – 61	60 – 65	64 – 70
cG Oberschenkelstrumpf mit Haftband extra weit	48 – 53	51 – 57	55 – 61	59 – 65	62 – 70	66 – 74	71 – 80



cG	44 – 53	46 – 57	50 – 61	53 – 65	56 – 70	60 – 74	64 – 80
cF	40 – 47	42 – 51	45 – 55	48 – 58	51 – 62	54 – 67	58 – 72
cE	30 – 37	33 – 40	36 – 42	38 – 44	40 – 46	42 – 48	44 – 50
cD	28 – 33	30 – 35	32 – 37	34 – 39	36 – 41	38 – 43	40 – 45
cC	29 – 34	31 – 36	33 – 38	35 – 40	37 – 42	39 – 44	41 – 46
cB ₁	24 – 26	26 – 28	28 – 31	30 – 33	32 – 35	34 – 37	36 – 39
cB	18 – 20	20 – 22	22 – 24	24 – 26	26 – 28	28 – 30	30 – 32
cY	28 – 30	29 – 32	31 – 34	33 – 36	34 – 38	36 – 40	37 – 41
cA	18 – 21	20 – 23	22 – 25	24 – 27	26 – 29	28 – 31	30 – 33

Länge	
/G normal	72 – 83
/G kurz	62 – 71
/D normal	38 – 43
/D kurz	33 – 37

Fussteil
Standard mit offener Spitze
Kleine Fussgrößen bis 23 cm: 36 – 38
Normale Fussgrößen ab 23 cm: 38 – 43