


## MASSTABELLE

mediven cotton®, mediven elegance®, mediven for men®,  
 mediven comfort®, mediven plus®, duomed basic®

Standard ohne Naht							
<b>cG Oberschenkelstrumpf mit Haftband</b>	43 – 48	45 – 52	49 – 56	53 – 60	56 – 64	60 – 68	64 – 72
<b>cG Oberschenkelstrumpf mit Haftband extra weit</b>	49 – 57	53 – 62	57 – 67	61 – 72	65 – 77	69 – 81	73 – 85
Umfänge	I	II	III	IV	V	VI	VII



<b>cG</b>	43 – 57	45 – 62	49 – 67	53 – 72	56 – 77	60 – 81	64 – 85
<b>cF</b>	39 – 52	41 – 56	44 – 60	47 – 65	50 – 69	53 – 73	57 – 78
<b>cE</b>	30 – 37	33 – 40	35 – 43	37 – 45	39 – 48	41 – 51	42 – 53
<b>cD</b>	27 – 33	29 – 36	32 – 39	34 – 42	36 – 45	38 – 48	40 – 50
<b>cC</b>	28 – 34	30 – 37	33 – 40	35 – 43	37 – 46	39 – 49	41 – 51
<b>cB<sub>1</sub></b>	23 – 27	24 – 29	26 – 32	29 – 35	31 – 37	33 – 39	35 – 41
<b>cB</b>	18 – 20	20 – 22	22 – 24	24 – 26	26 – 28	28 – 30	30 – 32
<b>cY</b>	26 – 31	28 – 33	29 – 35	31 – 37	32 – 38	33 – 40	34 – 42
<b>cA</b>	17 – 22	19 – 24	21 – 26	23 – 29	25 – 32	27 – 34	28 – 35

Länge	
<b>IG normal</b>	72 – 83
<b>IG kurz</b>	62 – 71
<b>ID normal</b>	39 – 44
<b>ID kurz</b>	34 – 38

Maxi-Leibteil AT mediven plus		
<b>Grösse I – III</b>	T = ab 95 cm	H = ab 110 cm
<b>Grösse IV – V</b>	T = ab 120 cm	H = ab 135 cm
<b>Grösse VI – VII</b>	T = ab 140 cm	H = ab 145 cm



Cosanum AG  
 Brandstrasse 28  
 8952 Schlieren  
 Happy Customer Service Ortho:  
 Tel. 043 433 66 77

## MASSTABELLE mediven forte®

Standard ohne Naht							
cG Oberschenkelstrumpf mit Haftband	44 – 47	46 – 50	50 – 54	53 – 58	56 – 61	60 – 65	64 – 70
cG Oberschenkelstrumpf mit Haftband extra weit	48 – 53	51 – 57	55 – 61	59 – 65	62 – 70	66 – 74	71 – 80

Umfänge	I	II	III	IV	V	VI	VII
 cG	cG 44 – 53	46 – 57	50 – 61	53 – 65	56 – 70	60 – 74	64 – 80
cF	cF 40 – 47	42 – 51	45 – 55	48 – 58	51 – 62	54 – 67	58 – 72
cE	cE 30 – 37	33 – 40	36 – 42	38 – 44	40 – 46	42 – 48	44 – 50
cD	cD 28 – 33	30 – 35	32 – 37	34 – 39	36 – 41	38 – 43	40 – 45
cC	cC 29 – 34	31 – 36	33 – 38	35 – 40	37 – 42	39 – 44	41 – 46
cB <sub>1</sub>	cB <sub>1</sub> 24 – 26	26 – 28	28 – 31	30 – 33	32 – 35	34 – 37	36 – 39
cB	cB 18 – 20	20 – 22	22 – 24	24 – 26	26 – 28	28 – 30	30 – 32
cY	cY 28 – 30	29 – 32	31 – 34	33 – 36	34 – 38	36 – 40	37 – 41
cA	cA 18 – 21	20 – 23	22 – 25	24 – 27	26 – 29	28 – 31	30 – 33

Länge	
IG normal	72 – 83
IG kurz	62 – 71
ID normal	38 – 43
ID kurz	33 – 37

Fussteil
Standard mit offener Spitze
Kleine Fussgrößen bis 23 cm = Schuhgrösse 36 – 38
Normale Fussgrößen ab 23 cm = Schuhgrösse 38 – 43